



BEYOND Polus Whitening Accelerator with Premier Treatment Kits + + + + +

Beyond Dental & Health
281.277.4352
www.beyonddent.com

Description

BEYOND Polus Whitening Accelerator with Premier Treatment Kits is an in-office whitening system designed to whiten within three, 10-minute sessions. **BEYOND Polus Whitening Accelerator** is a multi-functional, halogen-powered whitening lamp with *LightBridge* technology (combines halogen and LED light technologies), LED curing light and spot whitening device, and optional low-level laser therapy treatment device. It is indicated for dual-arch whitening; single-tooth whitening is also available. The built-in LED curing light has dual-wavelength output to cure any dental composite. The **BEYOND II Premier Treatment Kits** with *RENEW Touch-up* whitening gel are for single-patient use and are designed to improve oral health while whitening, with a reduced risk of sensitivity. **BEYOND II Advanced Formula Whitening Gel** contains 35% hydrogen peroxide and the *RENEW Touch-up* whitening gel contains 6% hydrogen peroxide. *RENEW* is designed to further increase the initial whitening results (one to two weeks after initial treatment) or to maintain the results (three months after initial treatment) with an in-office appointment. Special features of the light include three light-output settings (gentle, normal and intense) to meet individualized treatment needs; advanced filtration system designed to reduce harmful ultraviolet and heat output; *Touch User Interface (TUI)* panel for intuitive programming and treatment tailoring; motorized height adjustment; wireless remote control with a patient call button; and advanced temperature monitor and cooling system. There is a 30-day money-back guarantee and three-year manufacturer's warranty against factory defect.



Consultants' Comments

- ◆ "I have some patients who are very sensitive - having a gentle light output is great."
- ◆ "Three, 10-minute rounds - the system is quick to use."
- ◆ "This system produces the most consistent results. Every patient lightened at least three Vita shades."
- ◆ "A readout on the light or on the display would be nice."
- ◆ "Add a brush to the kit for removal of product near gum line."

Patients' Comments

- ◆ "It was fast."
- ◆ "That light seems pretty high tech."
- ◆ "My teeth look great."
- ◆ "I have bleached before, and this time my teeth were not sensitive."
- ◆ "I did not notice any scent to the products."

Ratings: Excellent + + + + + Very Good + + + + Good + + +

THE DENTAL ADVISOR • 3110 West Liberty, Ann Arbor, Michigan 48103 • 800.347.1330 • info@dentaladvisor.com

© 2011 Dental Consultants, Inc. No unauthorized duplication or reprints may be made. Inquiries concerning duplication may be directed to the publisher.

BEYOND Polus Whitening Accelerator includes the following:

- *BEYOND Polus Whitening Accelerator* lamp
- Polus LED curing light and spot-whitening device
- *BEYOND II Premier Single Patient Whitening Kit*
- Protective goggles
- Wireless remote control
- Replacement halogen bulb
- Online training (unlimited)
- Complete marketing package

BEYOND Premier Whitening Kits each contain the following:

- 4g, dual-barrel syringe of *BEYOND II 35% Whitening Gel* (packaged separately)
- 2 g syringe of *BEYOND BlueSeal Gingival Protection* w/one metal and one plastic tip
- 2.6g syringe of *RENEW 6% Touch-up Whitening Gel*
- 3g jar of polishing sand
- Face protection cloth
- 2mL jar of lip protection cream
- Two cotton rolls, cotton buds
- Suction tip for removing gel
- 3g jar of after-treatment fluoride
- 1.2-ounce tube of *BEYOND PearlWhite Toothpaste*
- Multilingual instruction manual
- Post-whitening care instruction card
- *VITA Bleachguide 3D-MASTER* paper guide
- Cheek retractor (packaged separately)
- Patient consent form (packaged separately)

BEYOND Polus Whitening Accelerator with Premier Treatment Kits was evaluated by two dental offices in 40 patients. This whitening system received a 98% clinical rating.

Suggested Retail Cost

- \$ 3,895.00/unit
- \$ 85.00/treatment kit

Product Features

BEYOND Polus Whitening Accelerator with Premier Treatment Kits was rated excellent for ease of use of the light and gel, completeness of the kit, and time of procedure. Consultants reported better than average whitening. **BEYOND Polus Whitening Accelerator** light is on a sturdy wheeled base and includes some features not found in any other whitening light. The three output settings helped tailor the treatment to the patient based on anticipated level of sensitivity. The articulation of the light head made it easy to position over the patients' teeth. The **Premier Treatment Kits** were reported to be the most complete kits consultants had encountered - all accessories required for the appointment were included. Most of the patients reported that they did not experience discomfort during the whitening process. Of those patients who did experience some discomfort, they rated it moderate. Average whitening was three to seven shades on the Vita guide with all of the patients reporting that their teeth were whiter. Sixty-seven percent of consultants reported **BEYOND Polus Whitening Accelerator with Premier Treatment Kits** was better than their current in-office tooth whitening system, whereas 33% reported it was equivalent. All consultants would switch to **BEYOND Polus Whitening Accelerator with Premier Treatment Kits** and recommend it to their colleagues.

Clinical Tips

- Keep *BEYOND II Advanced Formula Whitening Gel* refrigerated.
- Be sure to use the Face Protection Cloths to protect the lips from the light.
- Give the remote call button to the patient to call for assistance during the procedure.

Editors' note: *The manufacturer reports that the BEYOND BlueSeal Gingival Protection gel has an orange scent.* ■